


# Ag360™

## Antimicrobial silver coating



Experience the power of Ag360™, Precision Coating's proprietary, highly-adherent **silver** coating designed for superior performance. Clinically proven in medical implants, Ag360™ offers broad-spectrum **antimicrobial** protection with a low elution profile. Ag360™ demonstrates a significant reduction in **bacterial** growth and infections rates on applications such as vascular graft prostheses and dialysis catheters.

### Minimizes Bacterial Colonization



### Features & Benefits

- Broad-spectrum antimicrobial
- Highly adherent
- Long-lasting
- Minimal leaching
- Biocompatible

*K.S. Tweden, J.D. Cameron, A.J. Razzouk, R.W. Bianco, W.R. Holmberg, R.J. Bricault, J.E. Barry, and E.J. Tobin, "Silver Modification of Cardiovascular Textiles for antimicrobial Properties," ASAIO Journal, 43, M475-M481 (1997).*

### Applications

- Vascular grafts
- Active Implants
- Orthopedic implants
- Fixation devices
- Vascular access
- Dental devices

### Substrates

- PET
- Silicone
- PEEK
- Stainless steel
- Cobalt-chromium
- Titanium

Ag360™ is your solution for unparalleled **performance** and **infection control**.